



## Nipple Flow Rate Guide for Medical Professionals

Dr. Brown's® Nipples, when used in conjunction with the Dr. Brown's® Zero-Resistance™ Bottle Systems, provide reliable and consistent nipple flow rates to custom-fit an infant's feeding needs. Use this guide to help support clinical decisions.

Nipple Flow Rate Guide for Medical Professionals



Flow rate: 4.92 mL/min\*

### Ultra-Premie™ Nipple

- Excellent first choice for very premature or fragile infants.
- Appropriate for use when infant shows stress signals with the Premie Flow™ nipple.

Available in Narrow Only



Flow rate: 7.22 mL/min\*

### Premie Flow™ Nipple

- Frequently used for first oral feeding for premature infants or when infant shows stress signals while feeding with faster flow nipple level.
- Good choice for supplemental bottle feedings for infants feeding at the breast.

Available in Narrow and Wide-Neck



Flow rate: 9.88 mL/min\*

### Level T Nipple

- Excellent choice for infants transitioning from Premie Flow™ Nipple.
- Useful for late preterm or full-term infants.
- Another good choice for supplemental bottle feedings for infants feeding at the breast.

Available in Narrow Only



Flow rate: 13.31 mL/min\*

### Level 1 Nipple

- For more mature infant demonstrating consistent, efficient feeding skills.
- If considering for use with infant that is breastfeeding, please consult with a lactation consultant or feeding specialist.

Available in Narrow and Wide-Neck



**Internal Vent System is ESSENTIAL to maximize nipple performance.**

\*Pados, B. Milk flow rates from bottle nipples: What we know and why it matters. *Nursing for Women's Health*. 2021;25(3):229-235. <https://doi.org/10.1016/j.nwh.2021.03.006>



*Positive Feeding Experiences from the Start™*

**Zero-Resistance™ Bottle System Video**



**Dr. Brown's® Nipples**



**Narrow**



**Wide-Neck**

**Purchasing Information for Parents/Caregivers**



**Where to Buy - US**



**Where to Buy - Ca**

**medinfo@drbrownsmedical.com**