

# Essential Feeding Solutions

## Transition from Disposable to Reusable

Every NICU Baby deserves ESSENTIAL feeding tools for positive feeding experiences.



SB260-MED

Zero-Resistance™

### 60 cc AccuFeed Bottle System

- Begin with a slow flow nipple for initial safe feeding experiences.
- Graduations from 5 cc to 60 cc.
- Ultra-Premie™ option (SB266-MED) available.



SB411-MED

Zero-Resistance™

### 4 oz/ 120 mL Narrow Bottle System with Premie Flow™ Nipple

- The "transition-to-home" bottle for volumes over 60 cc.
- Level T (SB41606-MED) and Level 1 (SB413-MED) options available.



SN615-MED

### Narrow Ultra-Premie™ Nipple with Collar

- Slowest flow rate of the Dr. Brown's® Medical Nipple line.
- Excellent first choice for very premature or medically fragile infants.
- Appropriate for use when infant shows stress signals with Premie Flow™ nipple.



SN715-MED

### Narrow Level T Nipple with Collar

- Excellent choice for infants transitioning from Premie Flow™ Nipple.
- Useful for late preterm or full-term infants.



AC076-MED

### AccuBrush™

- Bristles designed to clean interior of AccuFeed bottle.
- Chevron cleaning tip at base of handle prevents damage to nipple during cleaning.



AC059-MED

### Non-Metal Bottle Brush

- Bristles designed to clean interior of 4 oz/120 mL or 2 oz/20 mL bottle systems.
- Chevron cleaning tip at base of handle prevents damage to nipple during cleaning.



AC032-MED

### Microwave Steam Sanitizing Bags

- Sanitize at least once every 24 hours; 20 cycles per bag.



SR101-MED

### 4.4" Wire Cleaning Brush

- Use to clean internal vent system.
- Not for use to clean nipple.
- Replace every 24 hours.

**Time to Make the Switch to REUSABLE**

*Positive Feeding Experiences from the Start.™*